

Different approaches to humidity generators, using concepts with carriergas-free evaporators cold trap and mixing design

Braunschweig 16.11.2017

Dr.-Ing.Still

ADROP Feuchtemeßtechnik GmbH
Kurgartenstraße 59
D - 90762 Fürth

Content:

- introduction company
- humidity generation with carriergas-free evaporators, examples
- humidity generation with cold trap, example
- humidity generation with gas mixing design, example

aDROP GmbH:


- start in 1997 as spin-off from the Department of Fluid Mechanics Erlangen (Prof. Durst).
- change of the legal form 1998 into 'aDROP Feuchtemeßtechnik GmbH'.
- key activities:
 - direct evaporators without carrier-gas and low dead volume [0,1 g/h ... 10 kg/h]
 - humidity generators (mainly custom design)
 - custom assemblies for research and production
 - contract research
- CAD 3D design and in-house manufacture, mostly in stainless steel
- microcontroller and electronic development (AVR 8bit, ARM 32bit)
- staff 8 coworkers + freelancers (CNC, LabView,..)

A: humidity generator with carriergas-free evaporator

direct evaporator:

- The fluid is dosed into the evaporator, which converts it into pure gas (overheated steam), without any additional (carrier-) gas.
- evaporation takes place within a porous stainless steel matrix. A temperature control supplies the necessary enthalpy. The interaction of pressure drop, temperature control and heatexchange provides a wide operational range, going up to 1:50 from min to max flow.
- high accuracy dosing in the liquid phase (eg. Coriolis massflow- controller) is much easier than accurate dosing in the gas phase.
- a stand-alone direct evaporator allows a modular rig design: Unlike spray gas evaporation, gas flow, steam flow, pressure and temperature can be adjusted individually.

A: humidity generator with carriergas-free evaporator




A: humidity generator with carriergas-free evaporator

- evaporation of fuel, up to 480°C evaporation temperature, additional gas on demand
- bench design with external mass flow controllers





A: humidity generator with carriergas-free evaporator

- calibration set-up for industrial high-temperature psychrometers,
- with test chamber up to 350°C



- A: humidity generator with carriergas-free evaporator

test chamber



generator with test chamber


•A: humidity generator with carriergas-free evaporator

- twin test rig for fuel cell testing
- 2x DV2-evaporator with dosing
- sequential control and DAQ with LabView
- including the safety provisions for H₂ and CO



A: humidity generator with carriergas-free evaporator

- two parallel evaporators, with a steam flow range of 1:1000 (2,5 g/h ... 2.500 g/h H₂O)
- T_{out} bis 600°C




working ranges of direct evaporator (referring to H₂O)

	Dynamik	Bereich Min.	Bereich Max.	Druck stand.	Druck Max.	Temp. Stand.	Temp. Max.
Direkt-Verdampfer							
μSTEAM	1:20 (50)	0,1...2 g/h	0,5...10 g/h	6 barg	100+ barg	200°C	
DV2	1:40 (50)	1...50 g/h	2,5...100 g/h	6 barg	100+ barg	200°C	400°C
DV4	1:40 (50)	10...500 g/h	50...2.500 g/h	6 barg	100+ barg	200°C	450°C
DV1C	1:40 (50)		100...4.000 g/h	6 barg	50 barg	300°C	450°C
DV3	1:40 (50)	100...6.000 g/h	200...10.000 g/h	6 barg	30 barg	300°C	400°C
Direkt-Verdampfer							
aTHMOS	1:10	100...1.000 g/h	200...2.000 g/h	6 barg		200°C	
aTHMOS-RS	1:10	500...5.000 g/h	600...6.000 g/h	6 barg		200°C	
Spray-Verdampfer (zB. SiCl₄)							
SGV1	1:50+	10...500 g/h	50...2.500 g/h	6 barg		250°C	450°C

B: humidity generator using a cold trap

cold trap:

- a humidified gas is cooled down below dewpoint, condensating the surplus water vapour
- peltier-operated cooling and precise temperature/pressure measurement can lead to high accuracy dewpoints.
- dewpoint range can be expanded by dilution provisions




B: humidity generator using a cold trap


- provides a stable rh- and temp.-environment for proteine cristalls while x-ray diffraction meas.
- Tdp and t stable within 0,03°C, range 3...40°C Tdp without pressure of dilution measures



aCCU-Feuchtegenerator
Oben: Gasversorgung
Unten: Steuer- und
Leistungselektronik




Meßkopf / Wärmetauscher mit
Goniometeraufnahme für
Proteinmontage in der
Röntgenkristallographieanlage



B: humidity generator using a cold trap

- temperature controlled test chamber for the humidity generator
- used for testing of concrete elongation while ageing etc.



Kammerdeckel geschlossen, mit Sensordurchführungen und Temperierungsleitungen




Kammerdeckel aufgeschoben mit innenliegender Probenhalterung



Kammer mit außenliegenden Temperierungsrohren und Isolation

C: humidity generator using a mixing design

- humidity generation using a reference sensor
- less severe demands regarding the humidity generation due to the presence of a control sensor
- high accuracy requires high precision control sensors



C: humidity generator using a mixing design

- aGEPRO with rated flow 0,2 ... 1,5 l/min
- T [5 ... 70°C], rh [5 ... 95%rF], accuracy control sensor ca. 1% rh, 0,2°C
- control mode: set T+rh or Tdp (dewpoint)



C: humidity generator using a mixing design

- aGEPRO with dewpoint mirror control for production
- small amounts of gases like NO_2 can be introduced to check cross sensitivity




- C: humidity generator using a mixing design



something different...

- generation of steam pulses with valve and nozzle
- using a high speed valve f_{rep} ca. 25 Hz
- using PWM modulation
- defined pulses below 1 sec
- continuous ca. 1,5 kg/h
- short time up to 5 kg/h

aDROP



ADROP Feuchtemeßtechnik GmbH,
Kurgartenstraße 59, D-90762 Fürth